

User Guide

Winchester Couch Drawer – COU92/DRW

DOHERTY

Thank you for purchasing a drawer for your Winchester Couch. Before installing your new drawer(s) please read these instructions carefully and ensure you fully understand them.

Kit of parts:




- x1 Drawer
- x4 Washers
- x4 Screws
- x1 Allen key

Specifications:

- | | | | |
|------------------|-----------------|------------------------|-----------------|
| • Width: | 567mm (22.3") | • Net weight: | 11.5kg (25.4lb) |
| • Height: | 132mm (5.2") | • Max load per drawer: | 10kg (22lb) |
| • Overall Depth: | 687.5mm (27.1") | | |





The Winchester couch allows the drawers to be positioned in multiple orientations. If the drawers are to be stacked underneath each other no more than a maximum of 3 drawers may be stacked together, as shown below.



Assembly:



1. Carefully remove all packaging and dispose of it in a sensible and environmental fashion.
2. Turn your Winchester Couch over and carefully place it upside down, ensuring the couch top is protected against dirt or sharp objects. **Sidhil Ltd. recommend that the couch is turned over with a minimum of two able bodied people.**
3. If you wish to attach a drawer at the foot end and the couch has a paper roll holder attached ensure the paper roll is facing outwards, if it is not remove the holder and turn it around.
4. Remove the inner drawer from the outer and place the outer on the couch frame, ensuring the large cut out section is facing upwards. Consider which side you want to access the drawers from when placing the outer drawer on the couch frame.
5. The drawer can be positioned at the head, middle or foot end. Select the desired position and align the slots in the outer drawer with the threaded inserts in the frame.
6. With the aid of the allen key fasten the outer drawer in place using the supplied washers and screws. The outer drawer is secured in 4 corners.

Correct orientation of paper roll holder if drawer to be placed at foot end



(Move on to point 9 if drawers are not to be stacked)

7. If you have purchased more than 1 drawer and you wish to stack them place the next outer drawer on top of the drawer below. Align the slots in the drawer with the threaded inserts in the drawer outer below.
8. As before, with the aid of the allen key fasten the outer drawer in place using the supplied washers and screws. **No more than a maximum of 3 drawers may be stacked together.**
9. Repeat the above processes until all required drawers are fastened in place. **Ensure all fasteners are tight.**
10. Turn the couch over and carefully position it back onto its feet. **Sidhil Ltd. recommend that the couch is turned over with a minimum of two able bodied people.**
11. Place all inner drawers back in the outer frames.
12. Your couch drawers are now assembled.



When using the drawers do not open more than 1 drawer at a time

Cleaning:

Routine cleaning for infection control must be carried out in accordance with your local infection control policy or regulatory body.

In the absence of this, all surfaces to be wiped down with a cleaning cloth, moistened with a mild detergent and diluted in hot water (60°C). Extra care should be taken around areas where excess dirt or dust may gather. Rinse down with clean water and dry off with a paper towel.

Always ensure the drawers (and couch) are dry before re-using.

Sidhil Ltd. does not recommend the use of concentrated bleach as this can mark painted surfaces.

Inspection:

The Winchester Drawers require very little maintenance, however if there are signs of damage or the drawers are not performing as they should, please contact your supplier for assistance.

Inspection points:

- Check there are no loose or missing parts or fixings
- Check fixings periodically (Sidhil Ltd. recommend every 6 months)

Please refer to your Winchester Couch user guide for warranty information in regards to this product.



Sidhil Limited Sidhil Business Park Holmfield Halifax West Yorkshire HX2 9TN

Tel: +44 (0)1422 233030 Fax: +44 (0)1422 233040 Email: info@sidhil.com
Sales Tel: +44 (0)1422 233030 Sales Fax: +44 (0)1422 233040

A Member of the Siddall & Hilton Ltd Group of Companies



This document was created with Win2PDF available at <http://www.daneprairie.com>.
The unregistered version of Win2PDF is for evaluation or non-commercial use only.