

MINIMUM STANDALONE/CLIENT SPECIFICATIONS

Operating Systems	Microsoft® Windows® 7 Professional (32-bit* or 64-bit) Microsoft Windows 10 Professional (64-bit) Note: If installing ELI™ Link with Diagnostic Cardiology Suite, system
Compatible Viewing Resolutions and Font Scaling	Resolutions
	Font scaling
	1024x768, 1366x768, 1280x800 100%
	1680x1050, 1920x1080, 1920x1200 100%, 125%, 150%
	2736x1824 100%, 125%, 150%, 200%
Processor	Performance equivalent or better than an Intel® Core™ i3
RAM	Minimum 4 GB
Hard Drive Capacity	Minimum 250 GB (If client: minimum 160 GB)
Archive	Network or external USB drive
Software Installation	USB
Network	100 Mbps connection or better
USB Ports	Requires at least two USB 2.0 ports (WAM/AM12 modules, external USB drive)
Input Devices	Standard keyboard, mouse, ECG acquisition input device (e.g. WAM, AM12 modules)
Printing Device	HP M501dn Windows Printer with HPUPD PCL5 driver or equivalent

MINIMUM SERVER REQUIREMENTS

Instrument Type	Performance equivalent to an Intel Xeon® class, quad-core with hyperthreading
Resolution	1024x768 or better
RAM	Minimum 4 GB
Operating System	Microsoft Windows 2012 R2 (64-bit) Microsoft Windows Server 2016
System Disk	100 GB for OS and product installation (RAID recommended for data redundancy)
Data Disks	550 GB hard drive space available HD controller with 128 MB read/write cache (RAID recommended for data redundancy)